


## Portable Raman for Quantification of Methanol in Contaminated Spirits

### Introduction

Over the past several years, an alarming trend has become evident, highlighting serious issues related to harmful illegally brewed alcohol globally. Home-distilled alcohol brewed from industrial solvents (i.e wood alcohol) and presented as an alcoholic beverage often contains methanol, which causes blindness and can lead to death when ingested. This issue has been reported with fatal consequences on several continents. (1-3) It came to a breaking point in September of 2012 when the Czech Republic officially banned the sale of hard liquor after 20 people died from the consumption of spirits with dangerous levels of methanol, which can occur when home brewing (2). After an exhaustive study of different screening tools, the Czech Republic turned to the use of portable Raman spectroscopy as the screening tool of choice for the identification and quantification of methanol in contaminated spirits. In this application note we will discuss the reasons why portable Raman spectroscopy is the ideal choice for this application and we will show a real world example of methanol-laced rum.

### Key Advantages


Raman is well recognized as an analytical tool for the discrimination of similar molecules such as ethanol ( $C_2H_5OH$ ) and methanol ( $CH_3OH$ ) as shown in Figure 1. But what makes Raman ideally suited for this particular application over comparative technologies such as FTIR is the ability to measure through optically transparent containers and its insensitivity to interference from water. These two key properties have allowed screeners to accurately test for the presence of methanol down to ~1% by volume in the field without the need to open the bottles being tested.



**Figure 1 Raman spectra of reagent grade ethanol (red) and methanol (blue).**

### Experiment and Results

To show the efficacy of this technique, we performed an in-house study measuring commercially available coconut rum which was diluted with methanol in concentrations between 0.33% and 5.36%. Similar to the method currently being used in Europe, we utilized the portable i-Raman® Plus because of its high sensitivity, high resolution, and fiber-optic probe sampling, to measure the Raman spectra of the mixtures as shown in Figure 2a. The spectra were then analyzed in BWIQ® chemometric software. A PLS regression model was developed on the data that was baseline corrected using airPLS. The two factor model developed over the range from 920  $\text{cm}^{-1}$  and 1580 $\text{cm}^{-1}$  gave the calibration curve shown in Figure 2b, which has a RMSECV of 0.17.



**Figure 2 Raman spectra of methanol laced rum with varying concentrations of methanol (a) concentration curve generated from PLS regression model (b).**

### Conclusion

These results verify the validity that portable Raman spectrometers can be used for not only the identification of adulterants, but also for providing quantitative results. This technique can be expanded to other types of adulteration including diethylene glycol- laced glycerin (4) and even water dilution of spirits.



## **References**

- (1) D. W. Lachenmeier, K. Schoeberl, F. Kanteres, T. Kuballa, E. M. Sohnius, J., *Addiction*, 106 (Suppl.1):20-30 (2011).
- (2) D. Spritzer and D. Bilefsky, *New York Times* 17 September 2012.
- (3) B. Collins, "Methanol poisoning: the dangers of distilling spirits at home", ABC Kimberley June 13, 2013, <http://www.abc.net.au/local/audio/2013/06/13/3781104.htm>
- (4) C. M. Gryniewicz-Ruzicka, S. Arzhantsev, L. N. Pelster, B. J. Westenberger, L. F. Buhse, and J.F. Kauffman, *Applied Spectroscopy*, 65(3):334-41 (2011).